CERTEST *Campylobacter*
Rapid One step *Campylobacter* spp. Card test

INTENDED USE

The CerTest *Campylobacter* is a one step coloured chromatographic immunoassay for the qualitative detection of *Campylobacter* spp. in faeces.

INTRODUCTION

Campylobacteriosis is the disease caused by the presence of *Campylobacter* spp. The common routes of transmission are fecal-oral, person-to-person sexual contact, ingestion of contaminated food or water, and the eating of raw meat. The onset of disease symptoms usually occurs two to five days after infection, but can range from one to ten days.

There are 16 species and six subspecies assigned to the genus *Campylobacter*, of which the most frequently reported in human disease are *C. jejuni* (subspecies *jejuni*) and *C. coli* (99% *C. jejuni*). *C. lariidis* and *C. upsaliensis* are also regarded as primary pathogens, but are generally reported far less frequently in cases of human disease.

The most common clinical symptoms of *Campylobacter* infections include diarrhoea (frequently with blood in the faeces), abdominal pain, fever, headache, nausea, and/or vomiting. The symptoms typically last three to six days.

PRINCIPLE OF THE TEST

The CerTest *Campylobacter* is a qualitative immunochromatographic assay for the determination of *Campylobacter* in fecal samples. The membrane is pre-coated with monoclonal antibodies, on the test band region, against *Campylobacter* antigens. During testing, the sample is allowed to react with the coloured conjugate (anti-*Campylobacter* monoclonal antibodies-red polystyrene microspheres) which was pre-dried on the test strip. The mixture then moves upward on the membrane by capillary action. As the sample flows through the test membrane, the coloured particles migrate. In the case of a positive result the specific antibodies present on the membrane will capture the coloured conjugate. The mixture continues to move across the membrane to the immobilized antibody placed in the control band region, a green coloured band always appears. The presence of this green band serves as 1) verification that sufficient volume is added, 2) that proper flow is obtained and 3) as an internal control for the reagents.

STORAGE AND STABILITY

Store as packaged in the kit at 2-30°C. The test is stable through the expiration date printed on the sealed pouch. The test must remain in the sealed pouch until use. Do not freeze.

PRECAUTIONS


- For professional *in vitro* diagnostic use only.
- Do not use after expiration date.
- All the specimens should be considered potentially hazardous and handled in the same manner as an infectious agent.
- The tests should be discarded in a proper biohazard container after testing.

SPECIMEN COLLECTION AND PREPARATION

Stool samples should be collected in clean containers and the assay should be done right after collection. The samples can be stored in the refrigerator (2-4 °C) for 1-2 days prior to testing. For longer storage, the specimen must be kept frozen at -20°C. In this case, the sample will be totally thawed, and brought to room temperature before testing. Ensure only the amount needed is thawed because of freezing and defrosting cycles are not recommended. Homogenise stool sample as thoroughly as possible prior to preparation.

Specimen preparation (see illustration):

- (1) Take out the cap of the stool collection tube and use the stick to pick up enough quantity of sample (approx. 150 mg), introduce the stick in 3 different parts of the sample. Add the sample into the stool collection tube. If the stool sample was liquid take 150 µl using a pipette.
- (2) Close the tube with the diluent and stool sample. Shake the tube in order to assure good sample dispersion.

MATERIALS PROVIDED

- CerTest *Campylobacter* card tests
- Instructions for use
- Stool collection tubes with sample diluent

MATERIALS REQUIRED BUT NO PROVIDED

- Specimen collection container
- Disposable gloves
- Timer


TEST PROCEDURE

Allow the tests, stool samples and controls to reach to room temperature (15-30°C) prior to testing. Do not open pouches until ready to perform the assay.


1. Proceed to shake the stool collection tube in order to assure good sample dispersion. Cut the end of the cap (3).
2. Remove the CerTest *Campylobacter* device from its sealed bag just before using.
3. Use a separate stool collection tube and device for each sample or control. Dispense exactly 4 drops or 100 µl into the circular window marked with an arrow, avoiding to add solid particles with the liquid (4).

In case the tests did not run due to solid particles fallen into the round window, stir the sample added or dispense a drop of extraction buffer until seeing the liquid running through the reaction zone.

- 4.- Read the result at **10 minutes** (the coloured bands appear).

INTERPRETATION OF RESULTS (please refer to the illustration below)

NEGATIVE: Only one GREEN band appears across the central window in the site marked with the letter C (control line).

POSITIVE: In addition to the GREEN control band, a RED band (test line) also appears in the site marked with the letter T (result line). That probably would mean *Campylobacter* infection.

INVALID: A total absence of the control coloured band regardless the appearance or not of the result line. Insufficient specimen volume, incorrect procedural techniques or deterioration of the reagents are the most likely reasons for control line failure. Review the procedure and repeat the test with a new test. If the problem persists, discontinue using the test kit and contact you local distributor.

NOTES ON THE INTERPRETATION OF RESULTS

The intensity of the red coloured band in the result line region (T) will vary depending on the concentration of antigens in the specimen. However, neither the quantitative value, nor the rate of increase in antigens can be determined by this qualitative test.

QUALITY CONTROL

Internal procedural controls are included in the test. A green line appearing in the control region (C) is the internal procedural control. It confirms sufficient specimen volume and correct procedural technique. A clear background is an internal negative background control. If the test is working properly, the background in the result area should be clear and not interfere with the ability to read the result.

EXPECTED VALUES

Campylobacter spp. are bacteria that are a major cause of diarrhoeal illness in humans and are generally regarded as the most common bacterial cause of gastroenteritis worldwide. In developed and developing countries, they cause more cases of diarrhoea than, for example, foodborne *Salmonella* bacteria. In developing countries, *Campylobacter* infections in children under the age of two years are especially frequent, sometimes resulting in death. In almost all developed countries, the incidence of human *Campylobacter* infections has been steadily increasing for several years. The reasons for this are unknown.

LIMITATIONS

1. The test must be carried out within 2 hours of opening the sealed bag.
2. An excess of sample could cause wrong results (brown bands appear). Dilute the sample with the buffer and repeat the test.
3. Some stool samples can decrease the intensity of the control line.
4. Positive frozen and thawed samples could cause wrong results. Use fresh stool samples.
5. This test provides a presumptive diagnosis of *Campylobacteriosis*. A confirmed infection diagnosis should only be made by a physician after all clinical and laboratory findings have been evaluated must be based in the correlation of the results with further clinical observations.

PERFORMANCE

Sensitivity and specificity

It was performed an evaluation of CerTest *Campylobacter*. It was studied 35 stool samples and the results were confirmed by ImmunoCard STAT! CAMPY. CerTest *Campylobacter* showed >99% of sensitivity and >99% of specificity.

The use of a mouse monoclonal antibody in CerTest *Campylobacter* assures high degree of specificity for the detection of these bacteria. The antibodies used to elaborate this test recognise epitopes found in stool of patients, as well as in preparations from the bacteria cultures *in vitro*.

This preliminary values has to be taken with precaution until more evaluation data will be available.

Cross-reactivity


It was performed an evaluation to determine the cross reactivity of CerTest *Campylobacter*. There is not cross reactivity with common intestinal pathogens, other organisms and substances occasionally present in feces:

- *H. pylori*
- *Escherichia coli* O157
- *Listeria monocytogenes*
- *Salmonella*

REFERENCES

- Fernández, H. and Farace, M.I. "Manual de Procedimientos *Campylobacter*". INEI. 2003.
- Kawatsu, K. et al. "Development and Evaluation of Immunochromatographic Assay for Simple and Rapid Detection of *Campylobacter jejuni* and *Campylobacter coli* in Human Stool Specimens". *Journal of Clinical Microbiology* Apr. 2008 Vol 46, No. 4, p. 1226-1231.

SYMBOLS FOR IVD COMPONENTS AND REAGENTS

CerTest BIOTEC S.L.
C/María de Luna 11, nave 16
E-50018 Zaragoza (SPAIN)
www.certest.es

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