CERTEST CALPROTECTIN 50+200

ONE STEP CALPROTECTIN
50+200 COMBO CARD TEST

CERTEST BIOTEC S.L.

ENGLISH

INTENDED USE

The CerTest Calprotectin 50+200 one step combo card is a coloured chromatographic immunoassay for the semi-quantitative detection of human fecal calprotectin (hFCP) in stool samples, to the following cut-off values: A: 50µg/g feces and B: 200µg/g feces.

INTRODUCTION

Calprotectin is a neutrophil cytosolic protein with antimicrobial properties, which is present at increased concentration in stool samples during bowel inflammation (intestinal inflammation). The stability of the protein to degradation keeps it stable in feces for up to 7 days at room temperature, making it an ideal analyte.

Calprotectin is released by activation of leukocytes, giving increased levels in plasma, cerebral spinal fluid, synovial fluid, urine or stools as a consequence of disease in the relevant organ(s).

Calprotectin inhibits zinc-dependent enzyme systems, as a result kills microbes and induces apoptosis in normal and cancer cells. In the presence of calcium, calprotectin is a remarkably resistant to proteolytic degradation.

This is a non-invasive marker of intestinal inflammation (for example in Ulcerative Colitis (UC) and Crohn's Disease (CD)).

CerTest Calprotectin 50+200 combo card offers a highly sensitive, simple and non-invasive assay to estimate inflammatory disease, to monitor treatment response in these patients, to predict risk of relapse and to monitor small bowel transplantation or graft rejection.

PRINCIPLE OF THE TEST

The CerTest Calprotectin 50+200 combo card is based on the principle of a semi-quantitative immunochromatographic assay for the determination of human calprotectin in stool samples. Each strip consists of: a membrane pre-coated with antibodies on the test line (result region), against human calprotectin and colored conjugate (anti-human calprotectin antibodies-red microspheres) pad pre-dried on the end of the membrane. During testing, the sample reacts with the conjugate and then the mixture moves upward on the membrane by capillary action. As the sample flows through the test membrane, the colored particles migrate. In the case of a positive result the specific antibodies present on the membrane will capture the colored conjugate. A colored band will be visible, in A and/or B depends on the calprotectin content of the sample. These bands are used to interpret the result.

The mixture continues to move across the membrane to the immobilized antibody placed in the control band region; this green colored band always appears. The presence of this green band serves as 1) verification that adequate volume is added, 2) proper flow obtained and 3) an internal control for the reagents.

STORAGE AND STABILITY

Store as packaged in the sealed pouch at 2-30°C. The test is stable until through the expiration date marked on its sealed pouch. The test must remain in the sealed pouch until use. Do not freeze.

PRECAUTIONS


- For professional *in vitro* diagnostic use.
- Do not use after expiration date.
- All the specimens should be considered potentially hazardous and handled in the same manner as an infectious agent. A new test must be used for each sample to avoid contamination errors.
- The tests should be discarded in a proper biohazard container after testing.

SPECIMEN COLLECTION AND PREPARATION

Stool samples should be collected in clean containers. The samples can be stored at room temperature or in the refrigerator (2-8°C) for 7 days prior to testing. For longer storage, maximum 6 months, the specimen must be kept frozen at -20°C. In this case, the sample will be totally thawed and brought to room temperature before testing.

Specimen preparation (see illustration):

- (1) Unscrew the cap and use the stick by introducing four times into the fecal specimen to pick up the sample and add the stick with the sample into the stool collection tube.
- (2) Close the tube with the diluent and stool sample. Shake the tube in order to assure good sample dispersion. The stool collection vial with dilute sample can be stored for 7 days in the refrigerator (2-8°C) prior to testing.

MATERIALS

MATERIALS PROVIDED

- CerTest Calprotectin 50+200 combo card
- Instructions for use
- Stool collection tubes with diluent


MATERIALS REQUIRED BUT NO PROVIDED

- Specimen collection container
- Disposable gloves
- Timer


TEST PROCEDURE

Allow the test, stool samples and controls to reach room temperature (15-30°C) prior to testing. Do not open pouches until the performing of the assay.


1. Proceed to shake the stool collection tube in order to assure good sample dispersion. Cut the end of the cap (3).
2. Remove the CerTest Calprotectin 50+200 combo card device from its sealed bag just before using.
3. Use a separate stool collection tube and device for each sample or control. Dispense 4 drops or 100µL in the circular window marked with the letter A (4) and 4 drops or 100µL, using the same tube, in the circular window marked with the letter B (5).
4. Read the result at **10 minutes** (the coloured bands appear). Do not read the test result after more than 10 minutes.

INTERPRETATION OF THE RESULTS (please refer to the illustration below)

INVALID

	A (50µg/g)	B (200µg/g)	Interpretation of the results
1.	-	-	Calprotectin marker is not present in patient sample (<50µg/g), which might mean neither active gastrointestinal inflammation, nor risk of relapse (CD or UC relapse), nor graft or transplantation rejection.
	GREEN	GREEN	
2.	+	+	High quantity of calprotectin marker (>200µg/g) is present in patient sample, which might mean active gastrointestinal inflammation, risk of relapse in clinical remission, graft rejection or transplantation rejection.
	GREEN-RED	GREEN-RED	
3.	+	-	Calprotectin marker is present in patient sample (>50µg/g and <200µg/g), which might mean active gastrointestinal inflammation.
	GREEN-RED	GREEN	
4.	-	+	Invalid result. If B is positive, A has to be also positive.
	GREEN	GREEN-RED	
5.	Any other result	Any other result	Invalid result: Invalid result either A or B, we recommend repeating the assay using the same sample with another combo device.

VALID RESULTS

- Strip A and Strip B: NEGATIVE.** Two GREEN bands appear across both central windows (A and B) in the site marked with the letter C (control line). That might mean neither active gastrointestinal inflammation is present in the patient, nor relapse risk is present in patients with Crohn's Disease and Ulcerative Colitis in clinical remission, nor graft or transplantation rejection
- Strip A: Calprotectin 50 POSITIVE.** Besides the GREEN control band, a second RED one (calprotectin 50 test line) appears in the site marked with the letter T (result region). **Strip B: Calprotectin 200 POSITIVE.** Besides the GREEN control band, a second RED one (calprotectin 200 test line) appears in the site marked with the letter T (result region). That might mean either active gastrointestinal inflammatory pathology, risk of relapse in patient during monitor treatment response (patients in clinical remission of GI inflammatory pathology, for example: Crohn's Disease or Ulcerative Colitis) or rejection in patient after small bowel transplantation or graft.
- Strip A: Calprotectin 50 POSITIVE.** Besides the GREEN control band, a RED band (calprotectin 50 test line) also appears in the site marked with the letter T (result region). **Strip B: Calprotectin 200 NEGATIVE.** A GREEN control band appears across the central window in the site marked with the letter C (control line). That might mean an active gastrointestinal inflammation in patient.

INVALID RESULTS

- Strip A: Calprotectin 50 NEGATIVE.** A GREEN control band appears across the central window in the site marked with the letter C (control line). **Strip B: Calprotectin 200 POSITIVE.** Besides the GREEN control band, a second RED one (calprotectin 200 test line) appears in the site marked with the letter T (result region).
- Any other result is considered to be invalid. A total absence of the control coloured band (GREEN) regardless the appearance or not of the result line (RED). → Insufficient specimen volume, incorrect procedural techniques or deterioration of the reagents are likely the reasons for control line failure. Review the procedure and repeat the assay using a new combo card test. If the problem still persists, discontinue using the test kit and contact your local distributor.

NOTES ON THE INTERPRETATION OF RESULTS

The intensity of the red colored band in the result line region (T) will vary depending on the concentration of antigens present in the specimen. However, neither the quantitative value, nor the rate of increase in antigen can be determined by this qualitative test.

QUALITY CONTROL

Internal procedural controls are included in the test. A GREEN line appearing in the control region (C) is an internal control, which confirms sufficient specimen volume and correct procedural technique.

LIMITATIONS

- The test must be carried out within 2 hours since the opening of the sealed bag.
- Stool from patients with non-steroidal anti-inflammatory drug (NSAID) treatment might show positive result.
- Neonatal fecal calprotectin levels have been reported higher than those in normal children with a mean of 167µg/g.
- An excess of stool sample could result in wrong results (brown bands appear or absence of the control colored band).
- Patients' stool suffering from active inflammatory bowel disease, usually involve significant neutrophilic inflammation of the intestine, such as Crohn's disease and ulcerative colitis, would be positive for fecal calprotectin. CerTest Calprotectin 50+200 combo card could be sensitive for this diagnosis in patients with chronic diarrhea.

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6. Positive results confirm the presence of human calprotectin in fecal samples. Nevertheless, it can also be due to several causes such as Inflammatory Bowel Disease, colorectal cancer and some other enteropathies. A positive result should be followed up with additional diagnostic invasive procedures, colonoscopy and biopsy in order to confirm the diagnosis and establish the inflammation extent.
7. Negative results do not exclude inflammation. Some diseases such as celiac sprue and microscopic colitis polyps may mainly involve mononuclear inflammation.

EXPECTED VALUES

Some studies established equal or higher 50µg hFCP/g feces as cut-off value to detect adult patients with gastrointestinal inflammatory pathology, requiring additional invasive procedures to diagnosis.

The clinical course of inflammatory bowel disease is characterized by a succession of relapses and remissions. Some studies established equal or higher 200µg hFCP/g feces as cut-off value to allow predict clinical relapse of disease activity in patients with Ulcerative Colitis and in Crohn's Disease and detect some neonatal patients (higher levels than normal children) with gastrointestinal inflammatory pathology that will require to diagnosis additional diagnostic invasive procedures.

Fecal calprotectin level evolution in individuals showed that fecal calprotectin level can predict rejection days before histopathological diagnosis.

PERFORMANCE CHARACTERISTICS

Analytical sensitivity: A sample containing calprotectin at concentration equal to or higher than 200µg/g feces produces positive results in two strip-test (A and B) when using CerTest Calprotectin 50+200 combo card. A sample containing calprotectin at concentration equal or higher than 50µg/g feces produces positive result when using strip A (calprotectin 50) of CerTest Calprotectin 50+200 combo card.

Different calprotectin dilutions were tested directly in the extraction buffer or spiked in a negative stool sample in accordance with the kit instructions to determinate the detection limit of the test.

Analytical specificity (Cross reaction): The CerTest Calprotectin 50+200 combo card test is specific for human calprotectin, showing no cross-reaction with other fecal markers: hemoglobin and transferrin.

Clinical sensitivity and specificity

64 stool patients' specimens were collected from two different hospitals and evaluated using CerTest Calprotectin 50+200 combo card. The detection of human calprotectin under CerTest Calprotectin 50+200 combo card showed >94% sensitivity correlation compared to another commercial immunoassay (Calprest® Eurospital).

The detection of human calprotectin using CerTest Calprotectin 50+200 combo card showed 93% specificity correlation compared to another commercial immunoassay (Calprest® Eurospital).

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1. Angriman I. et al. Enzymes in feces: Useful markers of chronic inflammatory bowel disease. *Clinica Chimica Acta* 381 Feb 2007, p. 63-68
2. Quail, M.A. et al. Fecal Calprotectin Complements Routine Laboratory Investigations in Diagnosing Childhood Inflammatory Bowel Disease. *Inflamm Bowel Dis*, Vol 15 No 5; May 2009, p. 756-759
3. Gaya D.R., et al. Faecal calprotectin in the assessment of Crohn's disease activity. *Q J Med* 2005, Vol 98, May 2005, p. 435-441.
4. Langhorst, M.D. et al. Noninvasive Markers in the Assessment of Intestinal Inflammation in Inflammatory Bowel Diseases: Performance of Fecal Lactoferrin, Calprotectin and PMN-Elastase, CRP, and Clinical Indices. *Am. J. Gastroenterol.* 2008; Vol 103, p. 162-169.
5. Tibble JA., et al., "Surrogate markers of intestinal inflammation are predictive of relapse in patients with inflammatory bowel disease". *Gastroenterology*, Vol 119 (1), Jul 2000, p. 15-22.
6. Gisbert J.P. et al. "Utilidad de la calprotectina y lactoferrina fecal en la predicción de la recidiva de la enfermedad inflamatoria intestinal (EII)". *Gastroenterología y Hepatología*, Vol. 32, Issue 3, March 2009, p. 226-227.

SYMBOLS FOR IVD COMPONENTS AND REAGENTS/

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