CERTEST CRYPTO-GIARDIA CARD
One Step *Cryptosporidium parvum* and *Giardia lamblia* Card test

INTENDED USE

The CerTest *Crypto-Giardia* Card is a one step coloured chromatographic immunoassay for the qualitative detection of *Cryptosporidium parvum* and *Giardia lamblia* in stool samples.

INTRODUCTION

Giardia and *Cryptosporidium* are parasites that can be found in water. *Giardia* causes an intestinal illness called giardiasis. *Cryptosporidium* is responsible for a similar illness called cryptosporidiosis. These infections have become the most common causes of waterborne diseases (found in both drinking and recreational water) in humans.

Giardiasis is a diarrheal illness caused by a very small parasite, *Giardia intestinalis* (also known as *Giardia lamblia*). Once an animal or person is infected with *Giardia*, the parasite lives in the intestine and is passed in the stool. The parasite is protected by an outer shell and can survive outside the body and in the environment for a long time. The most common symptoms of giardiasis include: diarrhea, loose or watery stool, stomach cramps and upset stomach.

Cryptosporidium parvum is the major cause of persistent diarrhoea in developing countries. This parasite is recognised as a highly infectious enteric pathogen and infective stage is transmitted by the fecal-oral route. Symptoms of cryptosporidiosis include watery diarrhoea, stomach cramps, weight loss, nausea and sometimes fever.

PRINCIPLE OF THE TEST

CerTest *Crypto-Giardia* Card uses a new one-step immunochromatographic system with coloured latex particles. It is a simple-to-use test which only needs the dilution of the fecal sample with the supplied ready-to-use dilution buffer.

Specificity is ensured by using mouse monoclonal antibodies directed against specific membrane antigens of *Giardia lamblia*. These reagents are used conjugated with coloured latex particles and coated on a nitrocellulose membrane. Liquid sample and coloured latex particles conjugate both migrate by capillary action. Upon reaching the first specific anti-*Giardia* monoclonal reagent, any *Giardia* present in the sample will be blocked and an immunoreaction will appear as a blue coloured line (first result line). The sample continues to migrate up the test, on reaching the second specific anti-Crypto monoclonal reagent, any *Cryptosporidium* present in the sample will be blocked and another immunoreaction will appear as a red coloured line (second result line). The sample continues to migrate up the test, on reaching the third reagent, a green coloured line (control line) is formed. This green line indicates that the sample has migrated correctly and that the chromatography has developed without hindrance. The control line appears with positive and negative samples. If this line does not develop, the test is invalid. The presence of this green band serves as verification that sufficient volume is added, that proper flow is obtained and as an internal control for the reagents.

STORAGE AND STABILITY

Store as packaged in the sealed pouch at 2-30°C. The test is stable through the expiration date printed on the sealed pouch. The test must remain in the sealed pouch until use. Do not freeze.

PRECAUTIONS


- For professional *in vitro* diagnostic use only.
- Do not use after expiration date.
- All the specimens should be considered potentially hazardous and handled in the same manner as an infectious agent.
- The tests should be discarded in a proper biohazard container after testing.

SPECIMEN COLLECTION AND PREPARATION

Stool samples should be collected in clean containers and the assay should be done right after collection. The samples can be stored in the refrigerator (2-4 °C) for 1-2 days prior to testing. For longer storage, maximum 1 year, the specimen must be kept frozen at -20°C. In this case, the sample will be totally thawed, and brought to room temperature before testing.

Specimen preparation (see illustration):

- (1) Unscrew the tap and use the stick to pick up a little sample, if the stool sample was liquid take 100 µL using a pipette, and add the sample into the stool collection tube.
- (2) Close the tube with the diluent and stool sample. Shake the tube in order to assure good sample dispersion.

MATERIALS PROVIDED

- Device tests
- Instructions for use
- Stool collection tubes with sample diluent


MATERIALS REQUIRED BUT NO PROVIDED

- Specimen collection container
- Disposable gloves
- Timer


TEST PROCEDURE

Allow the tests, stool samples and controls to reach to room temperature (15-30°C) prior to testing. Do not open pouches until ready to perform the assay.

1. Proceed to shake the stool collection tube in order to assure good sample dispersion. Cut the end of the top (3).
2. Remove the CerTest *Crypto-Giardia* Card device from its sealed bag just before using.
3. Use a separate stool collection tube and device for each sample or control. Dispense exactly 5 drops or 150 µL into the circular window marked with an arrow, avoiding to add solid particles with the liquid (4).
- 4.- Read the result at **10 minutes** (the coloured bands appear).

INTERPRETATION OF RESULTS (please refer to the illustration below)

NEGATIVE: only one GREEN band appears across the central window in the site marked with the letter C (control line).

CRYPTO POSITIVE: in addition to the GREEN control band, a RED band (Crypto test line) also appears in the site marked with the letter T (results lines).

GIARDIA POSITIVE: in addition to the GREEN control band, a BLUE band (Giardia test line) also appears in the site marked with the letter T (results lines).

CRYPTO-GIARDIA POSITIVE: All the lines above described (a GREEN control band in the control region, a RED band and a BLUE band in the result region) could appear at the same time during the test performance due to a simultaneous infection of Crypto and Giardia.

INVALID: A total absence of the control coloured band (GREEN) regardless the appearance or not of the results lines (RED/BLUE). Insufficient specimen volume, incorrect procedural techniques or deterioration of the reagents are the most likely reasons for control line failure. Review the procedure and repeat the test with a new test. If the problem persists, discontinue using the test kit and contact you local distributor.

NOTES ON THE INTERPRETATION OF RESULTS

The intensity of the red/blue coloured band in the result line region (T) will vary depending on the concentration of antigen present in the specimen. However, neither the quantitative value, nor the rate of increase in antigen can be determined by this qualitative test.

QUALITY CONTROL

Internal procedural controls are included in the test. A green line appearing in the control region (C) is the internal procedural control. It confirms sufficient specimen volume and correct procedural technique.

LIMITATIONS

1. The test must be carried out within 2 hours of opening the sealed bag.
2. An excess of stool sample could cause wrong results (brown bands appear). Dilute the sample with the buffer and repeat the test.
3. After one week of infection, the number of parasites in feces is decreasing, making the sample less reactive. Stool samples should be collected within one week of the onset of symptoms.
4. This test provides a presumptive diagnosis for Cryptosporidiosis and/or Giardiasis. A confirmed infection diagnosis should only be made by a physician after all clinical and laboratory findings have been evaluated.

EXPECTED VALUES

Negative results are expected in healthy patients. This test is intended to be used for the diagnosis of Cryptosporidiosis and Giardiasis.

PERFORMANCE

Sensitivity

The evaluation was conducted comparing the results obtained using the CerTest *Crypto-Giardia* Test to another commercial available ELISA *Crypto* and *Giardia* assays. The detection of *C. parvum* showed 95% of concordance in sensitivity. The detection of *G. lamblia* showed 95% of concordance in sensitivity.

Specificity

The use of mouse monoclonal antibodies in the elaboration of CerTest *Crypto-Giardia* Test assures its high degree of specificity for antigens of *C. parvum* and *G. lamblia*. This test shows >98% of specificity compared with commercial ELISA assays.

REFERENCES

1. MARSHALL, M.M., et al., "Waterborne Protozoan Pathogens", *Clinical Microbiology Review*, Jan. 1997, pp 67-85
2. DYLAN R. PILLAI and KEVIN C. KAIN, "Immunochromatographic Strip-Based Detection of *Entamoeba histolytica*-*E. dispar* and *Giardia lamblia* Coproantigen". *Journal of Clinical Microbiology*, Sept. 1999, Vol. 37, No 9, p. 3017-3019.
3. LYNNE S. GARCIA et al., "Commercial Assay for Detection of *Giardia lamblia* and *Cryptosporidium parvum* Antigens in Human Fecal Specimens by Rapid Solid-Phase Qualitative Immunochromatography", *Journal of Clinical Microbiology*, Jan. 2003, Vol. 41, No. 1, p. 209-212.

SYMBOLS FOR IVD COMPONENTS AND REAGENTS

	In vitro diagnostic device		Batch code
	Consult instructions for use		Catalogue number
	Keep dry		Contains sufficient for <n> tests
	Temperature limitation		Manufacturer
	Use by		Do not use if package damaged

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