INTENDED USE

The CerTest *Salmonella* is a one step coloured chromatographic immunoassay for the qualitative detection of *Salmonella* (*S. typhi*, *S. typhimurium*, *S. enteritidis*) in stool samples to detect a possible salmonellosis infection in humans.

INTRODUCTION

A total of 2501 different *Salmonella* serotypes have been identified up to 2004. *Salmonella enteritidis* and *Salmonella typhimurium*, the two most important serotypes for salmonellosis transmitted from animals to humans.

Salmonella infections are acquired by eating contaminated poultry, eggs or dairy products. The *Salmonella* spp. also spreads easily from raw or undercooked poultry to innocent vegetables, fruit or other foods via contaminated hands, knives, countertops or cutting boards. Due to the ability of *Salmonella* to multiply in a wide variety of foods, it is important to be able to isolate the organisms even when present in very small numbers in the faeces.

The symptoms of *Salmonella* infection are abdominal pain, diarrhoea, mild fever, chills, headache, nausea and vomiting, developing 12–72 hours (but occasionally as long as 7 days) after infection. The discomfort generally lasts a few days. It can be dangerous for the elderly, infants and the immunocompromised, who may become extremely ill. *Salmonella* is also one of the leading predictors for reactive arthritis, a painful, chronic and potentially debilitating condition that causes joint inflammation.

CerTest *Salmonella* provides a rapid detection of *Salmonella* spp. directly from the stool samples.

PRINCIPLE OF THE TEST

The CerTest *Salmonella* is a qualitative immunochromatographic assay for the determination of *Salmonella* (*S. typhi*, *S. typhimurium*, *S. enteritidis*) in stool samples. The membrane is pre-coated with antibodies, on the test band region, to recognize this antigen.

During testing, the sample is allowed to react with the coloured conjugate (anti-*Salmonella* antibodies-red microspheres) which was pre-dried on the test. The mixture then moves upward on the membrane by capillary action. As the sample flows through the test membrane, the coloured particles migrate. In the case of a positive result the specific antibodies present on the membrane will capture the coloured conjugate (result region). The mixture continues to move across the membrane to the immobilized antibody placed in the control band region, a GREEN coloured band always appears. The presence of this GREEN band serves as 1) verification that sufficient volume is added, 2) that proper flow is obtained and 3) as an internal control for the reagents.

STORAGE AND STABILITY

Store as packaged in the sealed pouch at 2-30°C. The test is stable through the expiration date printed on the sealed pouch. The test must remain in the sealed pouch until use. Do not freeze.

PRECAUTIONS


- For professional *in vitro* diagnostic use only.
- Do not use after expiration date.
- All the specimens should be considered potentially hazardous and handled in the same manner as an infectious agent.
- The tests should be discarded in a proper biohazard container after testing.

SPECIMEN COLLECTION AND PREPARATION

Stool samples should be collected in clean containers and the assay should be done right after collection. The samples can be stored in the refrigerator (2-4 °C) for 1-2 days prior to testing. For longer storage, the specimen must be kept frozen at -20°C. In this case, the sample will be totally thawed, and brought to room temperature before testing. Ensure only the amount needed is thawed because of freezing and defrosting cycles are not recommended. Homogenise stool sample as thoroughly as possible prior to preparation.

Specimen preparation (see illustration):

- (1) Unscrew the cap and use the stick to pick up the sample (approx. 150mg); introduce the stick in 3 different parts of the sample; if the stool sample was liquid, take 150 µL using a pipette, and add the sample into the stool collection tube.
- (2) Close the tube with the diluent and stool sample. Shake the tube in order to assure good sample dispersion.

Note: Because it is possible for the *Salmonella* content in the stool to be too small, a negative result cannot be considered to be the same as a negative antigen finding. It is recommended to carry out enrichment in Rappaport-Vassiliadis medium.

MATERIALS PROVIDED

- CerTest *Salmonella* card tests
- Instructions for use
- Stool collection tubes with sample diluent


MATERIALS REQUIRED BUT NO PROVIDED

- Specimen collection container
- Disposable gloves
- Timer


TEST PROCEDURE

Allow the tests, stool or enrichment samples and controls to reach to room temperature (15-30°C) prior to testing. Do not open pouches until ready to perform the assay.

1. Proceed to shake the stool collection tube in order to assure good sample dispersion. Cut the end of the cap (3).
2. Remove the CerTest *Salmonella* device from its sealed bag just before using.
3. Use a separate stool collection tube and device for each sample or control. Dispense 4 drops or 100 µL into the circular window marked with an arrow (4).
- 4.- Read the result at **10 minutes** (the coloured bands appear).

INTERPRETATION OF RESULTS (please refer to the illustration below)

- NEGATIVE:** Only one GREEN band appears across the central window in the site marked with the letter C (control line).
- POSITIVE:** In addition to the GREEN control band, a distinguishable RED band also appears in the site marked with the letter T (result line).
- INVALID:** A total absence of the control coloured band (GREEN) regardless of the appearance or not of the result line (RED). Insufficient specimen volume or incorrect procedural techniques are the most likely reasons for control line failure. Review the procedure and repeat the tests with a new test. If the problem persists, discontinue using the test kit and contact your local distributor.

NOTES ON THE INTERPRETATION OF RESULTS

The intensity of the red coloured band in the result line region (T) will vary depending on the concentration of antigen present in the specimen. However, neither the quantitative value, nor the rate of increase in antigen can be determined by this qualitative test.

QUALITY CONTROL

A green line appearing in the control region (C) is the internal procedural control. It confirms sufficient specimen volume and correct procedural technique.

EXPECTED VALUES

Salmonellosis constitutes a major public health burden and represents a significant cost to society in many countries. Although outbreaks usually attract media attention, studies indicate that more than 80% of all salmonellosis cases occur individually rather than as outbreaks. In the United States of America, an estimated 1.4 million non-typhoidal *Salmonella* infections, resulting in 168 000 visits to physicians, 15 000 hospitalizations and 580 deaths annually.

LIMITATIONS

1. The test must be carried out within 2 hours of opening the sealed bag.
2. An excess of sample could cause wrong results (brown bands appear). Dilute the sample with the buffer and repeat the test.
3. Some stool samples can decrease the intensity of the control green line.
4. Positive frozen and thawed samples could cause wrong results. Use fresh stool samples.
5. A negative result is not meaningful because it is possible the *Salmonella* content in the stool sample to be too small. A *Salmonella* determination should be carried out on a sample from an enrichment culture.
6. This test provides a presumptive diagnosis of *Salmonella* infections. A confirmed infection diagnosis should only be made by a physician after all clinical and laboratory findings have been evaluated must be based in the correlation of the results with further clinical observations.

PERFORMANCE

Detection limit
The detection limit for the different serotypes is: *S. enteritidis* 1x10⁴ bacteria/mL, *S. typhimurium* 1x10⁴ bacteria/mL and *S. typhi*: 1x10⁷ bacteria/mL.

Sensitivity and specificity
It was performed an evaluation using CerTest *Salmonella*. It was studied 40 stool samples and the results were confirmed by Singlepath@Salmonella. CerTest *Salmonella* test showed >99% of sensitivity and >97% of specificity. The antibodies used to elaborate this test recognise *Salmonella* epitopes found in stool of patients, as well as in preparations from the bacteria cultures *in vitro*. This preliminary values has to be taken with precaution until more evaluation data will be available.

Cross-reactivity

- It was performed an evaluation to determine the cross reactivity of CerTest *Salmonella*. There is not cross reactivity with common intestinal pathogens, other organisms and substances occasionally present in feces:
- *H. pylori*
 - *Escherichia coli* O157
 - *Listeria monocytogenes*
 - *Campylobacter*

REFERENCES

1. ABUBAKAR, I. "A systematic review of the clinical, public health and cost-effectiveness of rapid diagnostic tests for the detection and identification of bacterial intestinal pathogens in faeces and food". Health Technology Assessment 2007; Vol. 11: No. 36.

SYMBOLS FOR IVD COMPONENTS AND REAGENTS

	<i>In vitro</i> diagnostic device		Batch code
	Consult instructions for use		Catalogue number
	Keep dry		Contains sufficient for <n> tests
	Temperature limitation		Manufacturer
	Use by		Sample diluent

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